

How to activate your Evince or Breeze license

Table of contents

- How to start Evince and Breeze license wizard
- Activate Evince or Breeze license
 - Option 1: Using your username and password from Prediktera User pages
 - Option 2: Using license text that was emailed to you by Prediktera
- Generate activation code and send to Prediktera (Not required for trial)

Note. New users please first register an account at Prediktera user pages:


<http://user.prediktera.se>

Note: Evince and Breeze requires Java to run. Please follow instructions here after installation:


http://www.prediktera.se/download/pdf/Prediktera_Installing_Java.pdf

How to start Evince and Breeze license wizard

In Evince. Click on *Help* and *Edit License*




In Breeze. Right click on *Shutdown* button and *Edit license*



Activate Evince or Breeze license


Using your username and password from Prediktera User pages

1. Start license wizard. Select *Activate license by login* and click *Next*



The screenshot shows a dialog box titled "License Wizard - Step 1 of 2". The main heading is "Welcome to the Evince and Breeze license wizard". Below this, it says "Please choose how you wish to activate your license." and includes a note for new users to register at <http://www.prediktera.se>. There are two sections: "Online" and "Offline". Under "Online", the option "Activate license by login" is selected with a radio button, and it specifies "Username or email and Password from Prediktera user pages". Under "Offline", there are two unselected options: "Activate license from text" (sent to you by email) and "Generate activation code". At the bottom, there is a prompt "Click 'Next' to continue..." and a set of buttons: a help icon, "< Back", "Next >" (highlighted with a blue dashed border), "Finish", and "Cancel".

2. Enter your username and password and click Finish




The screenshot shows a dialog box titled "License Wizard - Step 2 of 2". The main heading is "Activate license by login". It states "Before you can continue, you must login. Please enter your username or email and password and press 'Finish'". It includes the same note for new users and the URL <http://www.prediktera.se>. There are two input fields: "Username or email" and "Password". A "Proxy Settings" button is located below the password field. At the bottom, there is a set of buttons: a help icon, "< Back", "Next >", "Finish", and "Cancel".


3. The License wizard closes automatically. Restart Evince or Breeze to start using your activated installation.

Using license text that was emailed to you by Prediktera


1. When you receive the license from Prediktera. Start the license wizard.
2. Select *Activate license from text* and click *Next*



3. Find the email that was sent to you from Prediktera when registering. Copy the whole license text.



4. Click *Paste from clipboard*. Click *Finish* to complete the license activation.




5. The License wizard closes automatically. Restart Evince or Breeze to start using your activated installation.


Generate activation code and send to Prediktera

This step is only required for activation of a full license (Not for trial)

1. Start license wizard. Select *Offline/Generate activation code* and *Next*



2. Click on *Copy to clipboard*



3. Paste the whole activation text in an email and send to license@prediktera.se
4. After you receive an email from Prediktera with the license, follow the instructions above "Using license text that was emailed to you by Prediktera"